

Job Description	Treat a Contaminated Tray and Prepare Instruments for Sterilization Participant will assemble materials and equipment necessary, separate waste, sort instruments, place into ultrasonic, turn on and verbalize process, remove instruments, and prepare for autoclaving.
Administration Time	10 minutes
Equipment, Supplies, or Materials Required	<ul style="list-style-type: none">• Sink• Hand soap/Hand sanitizer• Personal protective equipment:<ul style="list-style-type: none">○ utility gloves○ mask○ gown○ safety glasses• Instruments for autoclave<ul style="list-style-type: none">○ basic setup (mouth mirror, explorer, cotton pliers)○ high speed handpiece• Ultrasonic and solution• Autoclave/Nyclave bags/Autoclave wrap• Autoclave tape• Pencil or pen• Paper towels• Waste container• Gauze (made to appear contaminated with blood)• Other disposable items• Biohazard bag